Photometer PF-3

Handbuch
Manual
Manuel

Kompaktphotometer PF-3
Compact photometer PF-3
Photomètre compact PF-3
Fotómetro compacto PF-3

MACHEREY-NAGEL

www.mn-net.com

de

Handbuch PF-3 4

en

Manual PF-3 17

fr

Manuel PF-3 29

es

Manual PF-3 42

Inhalt

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sales@mn-net.com · www.mn-net.com

Contents

1. Introduction	18
1.1. Technical data	18
2. Safety precautions	18
2.1. Chemical/biological safety	19
2.2. Mechanical safety	19
2.3. Inappropriate use	19
3. Set-up	19
3.1. Package content	19
3.2. Instrument environment / set-up	19
3.3. Interface / power supply	19
4. Operation	20
4.1. Turning on the instrument	20
4.2. User guidance	20
4.3. Icon explanation	21
5. Measuring	21
5.1. Test selection	21
5.2. Additional tests	21
5.3. <i>VISOCOLOR</i> [®] and <i>NANOCOLOR</i> [®] tests	21
5.4. Measurement	22
5.4.1. Measurement range bar	22
5.4.2. 0-Set function	23
6. Settings	23
6.1. Activate / deactivate tests	23
6.2. Date	23
6.3. Time	24
6.4. Reaction time	24
6.5. Auto off	24
6.6. Update	24
6.7. Calibration	24
6.8. Measurement range bar	25
6.9. 0-Set function	25
6.10. Calibration adjustment	25
7. Memory	26
7.1. Display memory	26
7.2. Export data	26
8. User-defined special methods	26
9. Cleaning the instrument	26
10. Inspection equipment monitoring	27
11. Troubleshooting	27
12. Service	27
12.1. Spare parts and accessories	27
12.2. Disposal	27
12.3. Warranty	27
12.4. General usage	28
12.5. Contact	28

1. Introduction

Thank you for purchasing a PF-3 compact photometer from MACHEREY-NAGEL. The PF-3 is a quick and powerful instrument for mobile use. With the instrument, you can determine our *VISOCOLOR*[®] and *NANOCOLOR*[®] tests. Which tests can be run on the instrument depends on the version.

1.1. Technical data

Type:	LED photometer with microprocessors control, self-test and auto calibration
Optics:	LEDs + interference filters Insensitive to external light for quick measurements without light shield
Wavelength:	Version dependent, 3 wavelengths possible
Wavelength accuracy:	± 2 nm, half width 10–12 nm
Number of filters:	3
Light source:	LED
Detector:	Silicium photo element
Measurement modes:	Pre-programmed tests (<i>NANOCOLOR</i> [®] tube tests and <i>VISOCOLOR</i> [®] <i>ECO</i> tests)
Cuvette slot:	Tubes 16 mm OD
Memory:	50 results
Display:	Backlit graphic display, 64 x 128 px
User guidance:	High quality key foil, all data at a glance: result and unit, date, time
Interface:	Mini-USB
Update:	Free of charge via internet / PC
Operation environment:	5–50 °C at 90 % rel. humidity
Power supply:	3 AA batteries, rechargeable batteries, USB interface; optional internal accu-pack
Housing:	Water proof according to IP 68 (1 h, 1 m)
Dimensions:	170 x 95 x 68 mm
Weight:	0.5 kg
Warranty:	2 years

 The instrument complies with the following directives:
- 2006/95/EG - Low-Voltage Directive
- 2004/108/EG - EMV Directive

2. Safety precautions

Please read the manual carefully and completely before setting up and using the instrument. Please adhere to all remarks. Not complying with this information may cause harm to the user or damage the instrument. In order not to compromise the instruments safety measures, it may never be installed or used in any other way than described in this manual.

2.1. Chemical/biological safety

When using the instrument, it may be necessary to handle hazardous chemicals or biologically harmful samples. Please adhere to the following guidelines when using the instrument:

- Please read all information concerning hazards and safety measures in the material safety data sheets and on the packages of the solutions you work with.
- All used up solutions need to be disposed of in agreement with national rules and regulations.
- The type of safety equipment needs to be in line with the concentration and hazard potential of the given substance used.

2.2. Mechanical safety

If the housing is damaged, you need to send the instrument in for repairs. Instrument functions and results cannot be guaranteed anymore.

In case the USB cable is damaged while the power is supplied via USB, you need to stop using the cable and exchange it immediately.

2.3. Inappropriate use

The PF-3 is only meant to be used as a photometric analysis device for MACHEREY-NAGEL test kits. In case of inappropriate use by untrained staff, MACHEREY-NAGEL declines any warranty claim.

3. Set-up

3.1. Package content

Open the package carefully. If parts are missing or damaged, please contact your distributor or our technical service (*see 12.5, page 28*).

Depending on the Version, the instrument may come in a cardboard box or case. The case may contain reagents.

3.2. Instrument environment/set-up


Only set-up the instrument at a suitable location. The operation temperature ranges from 5–50 °C, humidity can be up to 90%. Protect the instrument from direct sunlight and strong temperature deviations.

3.3. Interface/power supply

The instrument features a USB port for power supply and data export. The interface is beneath the battery compartment cover at the back of the instrument and can be opened with a screw driver.

USB interface

Battery compartment

Initially, the instrument is equipped with 3 AA batteries. You can also use rechargeable batteries or an optional accu-pack (REF 919 391), which can be charged via the USB interface. It is also possible to run the instrument via the USB port only.

4. Operation

4.1. Turning on the instrument

Press the button with the On / Off symbol once. The instrument powers up and the screen is illuminated. You can see instrument name and firmware version. During the start-up phase, the instrument also performs a self-check. Then, the instrument jumps directly into the measurement mode.

4.2. User guidance


The instrument is designed for easy and intuitive use. The user guidance is based entirely on icons, which can be activated via the four activity-buttons below the display.

The screen itself is subdivided into three areas. The upper status bar shows date and time as well as the battery life.

In the middle of the screen you can see the actual data such as results or menu icons.

The lowest bar contains the activity-icons. The icons always represent a concrete function, which can be accessed with the underlying activity-buttons. When leaving menus via , your entries are saved. If you leave the menu via , none of your selections are activated or saved.

4.3. Icon explanation

	/		Scroll left / right		Escape
			Measure zero		Jump one level back
			Measure sample		Change number
			Confirm / Selection		Activate / Deactivate

5. Measuring

5.1. Test selection

To start a measurement, you need to select the respective test first. You can do so using the  /  arrows. Normally, all *VISOCOLOR*[®] tests, which are enabled for your PF-3 version, are activated on the upper menu level as favorites. These tests can be activated or deactivated (see 6.1, page 23). All *NANOCOLOR*[®] tests are initially deactivated, but can be activated any time via the settings menu.

Some tests feature multiple submethods to show different units. These submethods are saved within the instrument as a normal test. Thus, you can choose exactly the right unit / submethod that is relevant to your measurements.

Note: We recommend to assign your favorite measurements to the upper menu level when you first use the instrument, so that you can measure exactly the tests you need for your work quickly and easily.

5.2. Additional tests

In many cases, the wavelengths built into each version of the instrument enable determining additional tests and parameters, which are not included in the initial instrument version. If available, you can download such tests free of charge via the internet and add them with a small software tool to the instrument. You can find the respective software and a list with available tests at www.mn-net.com/PF-3.

5.3. *VISOCOLOR*[®] and *NANOCOLOR*[®] tests

The PF-3 is designed to run both *VISOCOLOR*[®] and *NANOCOLOR*[®] test kits. For the *VISOCOLOR*[®] tests, we provide a manual with pictogram test instructions. The manual is included in the case versions of the instrument and available for download under www.mn-net.com. The instructions on how to run the *NANOCOLOR*[®] tests can be found in the respective instructions of the kits.

An additional difference between the two test systems is that *VISOCOLOR*[®] kits always require measuring a blank value, while for almost all *NANOCOLOR*[®] tests the blank value is pre-stored in the instrument, thus enabling direct measurements. You can


also save one zero measurement for all *VISOCOLOR*[®] tests with our 0-Set function, so that you need to perform a zero measurement only once when running *VISOCOLOR*[®] test kits (see 5.4.2, page 23).

Apart from different parameters, *VISOCOLOR*[®] and *NANOCOLOR*[®] tests always marked with the addition of “NANO” or “VISO”, respectively and a 4-digit number. The number is unique for each test and refers to the REF number. For *NANOCOLOR*[®] tests, the number always starts with a 0, while *VISOCOLOR*[®] tests are indicated with a 5.

5.4. Measurement

To measure a test, prepare the sample (and the blank value if necessary) according to the test instructions. Then, perform the following steps (if there is no blank value necessary, you can ignore step 2–4).

1. Select test
2. Insert blank tube into the instrument
3. Press 
4. Remove blank tube
5. Insert tube with sample and reagents
6. Press 
7. The result is displayed
8. Once you remove the tube, the result is saved in the instrument memory

If the result is outside the test’s measurement range, the instrument shows < or > and the lowest value of the measuring range. If a measurement value was determined with a method, whose calibration was adjusted, the result will be displayed with an exclamation mark right beside the measurement value. For deleting the calibration adjustment of a method, the factory-provided calibration has to be recovered (see 6.10, page 25).

5.4.1. Measurement range bar

In addition to the result, the instrument can optionally display a measurement range bar, to check where the result lays in terms of the test’s measuring range. By default, the function is not active,

but you can turn it on or off, respectively in the settings (see 6.8, page 25).

5.4.2. 0-Set function

Generally, all VISOCOLOR® test require a zero measurement. Hence, if you switch tests, you need to perform a new zero measurement. If you run a test sequence for different parameters from the same sample, you can set one zero measurement for all subsequent VISOCOLOR® tests, which means you only need to measure a zero once.

If the 0-Set function is activated and you have performed a zero measurement, you can select a different test, insert the measuring tube directly and run the measurement via . If you measure a new zero, it will overwrite the first one and be stored as the next superordinate zero value.

Once you turn-off the instrument, the 0-Set value is lost and you need to measure a new zero.

By default, the function is not active, but you can turn it on or off, respectively in the settings

Note: We recommend to only use the same zero measurement as a basis for different measurements if the tests are performed from the same sample, or if the samples are very similar.

If you are not sure if the similarity is adequate, we always recommend measuring a new zero for our VISOCOLOR® tests.

6. Settings

You can find the settings menu on the upper menu level using the  /  arrows. To enter the settings menu, scroll to the corresponding icon and press the activity button .

6.1. Activate/deactivate tests

Within this menu, you can decide which tests are selected as favorites in the upper menu level. To scroll through the tests, use  / . Press  to activate or deactivate a test. The symbol in the upper right corner indicates whether the test is active or not ( / .

6.2. Date

In the date menu, you can set the date. As soon as you select the menu, you see the currently set date. The first number is highlighted and can be adjusted via the  button. With the  button you can scroll through the individual figures of the date.

Once you have set the date, you can chose the date format from the following options via  / .

- dd.mm.yyyy
- mm/dd/yyyy
- yyyy-mm-dd

Press the  activity button to save your changes. If you leave the menu via , your choices are discarded.

6.3. Time

In the time menu, you can set the time. As soon as you select the menu, you see the currently set time. The first number is highlighted and can be adjusted via the  button. With the  button on the right side you can scroll through the individual figures of the time.

Once you have set the time, you can adjust the time format via  / . You can choose between 24 h and 12 h mode.

Press the  activity button to save your changes. If you leave the menu via , your choices are discarded.

6.4. Reaction time

In the reaction time menu, you can decide whether the instrument should display the reaction time for the test you selected and serve as a counter function. If the reaction time is deactivated, you can start your measurement right away. You can activate / deactivate the reaction time using the  button.

Note: If the reaction time is active, press  twice to start a measurement ahead of time.

6.5. Auto off

In the auto-off menu, you can set if and when the instrument turns off after a certain time automatically. Use  /  to choose between 4 different auto-off times (5, 10, 15, 20 min) and continuous mode (∞).

Press the  activity button to save your changes. If you leave the menu via , your choices are discarded.

6.6. Update

In case you want to update the instrument, you need activate the update in the update menu. If there are updates available, you can find them online on our homepage under www.mn-net.com/PF-3. More detailed instruction regarding updates can be found on the web page in the respective update instructions.

6.7. Calibration

In the calibration menu, you can calibrate the instrument. Normally, the PF-3 does not need to be calibrated. In case a calibration is necessary, press  to start the calibration process. A new display appears asking you to insert the calibration tube. Press  to start

the calibration. Once the calibration is complete, you can leave the menu via .

6.8. Measurement range bar

As described in the measurement section, the instrument can optionally display a measurement range bar to check where the result lies with regards to the measuring range of the test. In the measuring range bar menu, you can activate / deactivate the bar using the  button.

6.9. 0-Set function

As described in the measurement section, the instrument can optionally store one zero value for all following *VISOCOLOR*[®] tests. In the 0-set function menu, you can activate / deactivate the function using the  button.

6.10. Calibration adjustment

The PF-3 device is provided with the MACHEREY-NAGEL factory calibration and can be readily used without any need for user calibration. The calibration adjustment tool can be used to meet regulatory requirements. Therefore the factory-provided calibration of the *VISOCOLOR*[®] and *NANOCOLOR*[®] methods can be adjusted with a standard in the menu “calibration adjustment”. You can choose the test to be adapted via . Confirm with  and insert the blank into the cuvette slot and press the button . Afterwards insert the reacted standard and press . The displayed value can now be changed with the button . The button  can be used to change the digit to be edited. The standard calibration adjustment can only be done within the measuring range of the test. Therefore the entry of the value zero is not allowed. Please confirm your entry with the button . The new button  will appear in the display, indicating that the calibration was adjusted for the test. To restore the factory-provided calibration, you have to confirm the button . The button will disappear and the test is working with the MACHEREY-NAGEL provided calibration again. After calibration adjustment or its removal, you can leave the menu via the button . If the calibration adjustment was done for a method, the results measured with this method will be indicated with an exclamation mark in the measurement display as well as in the memory.

Advice: We recommend to use the factory-provided calibration and not to change this one. The calibration can be adjusted to meet regulatory requirements. For the calibration adjustment procedure we recommend to use a standard with a value located in the upper third of the measurement range to reduce the calibration error to a minimum.

7. Memory

The PF-3 features an internal memory for up to 50 measurements. In addition to the actual result, the instrument also saves the date and time of the measurement. If the memory is full, the instrument overwrites the oldest data set once a new measurement is performed.

7.1. Display memory

Within the display memory menu, you can find all saved results. You can scroll between the results using  / . In the upper right hand corner of the screen, you see the number of the entry and the total number of entries within the memory. If a value in the memory is indicated with an exclamation mark, the value was measured with a method, whose factory-calibration was adjusted by the user (see 6.10, page 25).

7.2. Export data

You can export data from the instrument's memory to a PC. This option is only displayed in the screen if the instrument is actually connected to a computer via the USB port. To do so, connect the instrument to the USB port of your computer (see 3.3, page 19).

To transfer data, MACHEREY-NAGEL offers a special, free of charge software for convenient data export. You can download the respective software at www.mn-net.com/PF-3. You can find all additional information in the software instructions, which are also available on our website.

8. User-defined special methods

Besides the pre-programmed *VISOCOLOR*[®] und *NANOCOLOR*[®] tests, up to ten additional user-defined special methods can be stored in the PF-3 device. For the creation and administration of the user-defined special methods you need to install the *NANOCOLOR*[®] Data export software, which is available free of charge on the MACHEREY-NAGEL homepage via www.mn-net.com/PF-3. For further information please refer to the software manual of the data export software, which can be downloaded free of charge from our website.

9. Cleaning the instrument

Clean the housing of the instrument with a damp cloth. Wipe away splashes as soon as possible. Clean and dry the cuvette slot with a soft cotton tissue. Generally, tubes and photometer always need to be kept clean. Clean the tubes with a lint-free cloth.

10. Inspection equipment monitoring

The photometric accuracy of the PF-3 can be tested with the colored NANOCHECK solutions (REF 925 701), provided by MACHEREY-NAGEL. The appropriate evaluation sheet can be downloaded on our homepage www.mn-net.com.

11. Troubleshooting

Error message	Reason	Solution
Display stays dark	• No power supply	• Check if batteries are inserted or if USB cable is connected
Instrument is turned on but doesn't react	• Instrument is frozen	• Turn off instrument, wait 10 s, turn on instrument again

12. Service

12.1. Spare parts and accessories

Description	REF
Robust case with foam inlay, for PF-3, accessories and VISOCOLOR® test kits	934 102
Calibration tube	916 908
USB cable for data export and power supply	919 390
Optional accu pack	919 391
Measurement tubes (20)	916 80
Manual	919 392
VISOCOLOR® test instructions for PF-3	934 001

12.2. Disposal

In compliance with local and national legal regulations (EU Directive 2002/96/CE), the MACHEREY-NAGEL company disposes old instruments free of charge.

Note: With effect from August 2005, disposal using public waste disposal facilities is no longer permitted. In the case of disposal, please contact your MACHEREY-NAGEL representative

12.3. Warranty

The guarantee period is 24 months. The original commercial invoice serves as the guarantee card and must be presented in case of a claim during the given period. The guarantee is not valid in cases of improper handling and the guarantee does not cover defects caused by the external power supply. The guarantee is limited to repairs of defective parts or, in our sole discretion, replacement with a brand new and defect-free instrument. Any guarantee claims made or alleged do not influence the total guarantee period of 24 months. The right of withdrawal does not exist. Further claims are excluded and

these especially include claims for compensation caused by direct or indirect damages. Additionally our current Terms and Conditions for the Sale of Goods as printed on price lists shall apply.

12.4. General usage

Only use the instrument as instructed to analyze MACHEREY-NAGEL tests. Please do not open the instrument or make any other unauthorized modifications. The instrument is a highly sensitive measuring device with high accuracy for the photometric evaluation of MACHEREY-NAGEL tests. Unauthorized modifications or opening the instrument in an inappropriate way can result in maladjustments of the instrument's optical components, or similar damage, which in turn can lead to wrong results.

Note: Due to the reasons stated above, MACHEREY-NAGEL declines any responsibility for the correct functioning of the instrument if it has been opened or used improperly and objects to any warranty claims made in such a case. MN's warranty and liability obligations as to product(s) or to any party on account thereof shall be void and of no effect if the product(s) or any part thereof have been subject to unauthorized change, abuse, misuse, accident or alteration, unauthorized application or installation not in compliance with the product(s) specifications as set forth in the product leaflets / instructions hereof ("unauthorized action") and such unauthorized action resulted in the non-conformance or defect of the product(s).

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
sales@mn-net.com · www.mn-net.com

MACHEREY-NAGEL

Your competent partner in analytical chemistry

Filtration

Rapid Tests

Water Analysis

Chromatography

Bioanalysis

Tradition and Modernity – 100 years experience

- Worldwide operating German company (founded 1911)
- Subsidiaries in France, Switzerland and USA
- Distributors in over 150 countries
- Development, production and sales of special products for water, environmental and food analysis, for biotechnology, chemical and pharmaceutical industry and medical diagnostics

**Investment in the future –
New building in Düren for Research & Development,
production and warehouse**

Service

- Technical support by our Customer Service Center
- Product-specific quality certificates and Material Safety Data Sheets
- Customized special products
- Easy ordering and short delivery times
- Large number of catalogs, flyer and product information – of course free of charge

Contact

Technical support and customer service

Filtration, Rapid Tests, Water Analysis: +49 24 21 969-161
+49 24 21 969-138
+49 24 21 969-190
+49 24 21 969-174
+49 24 21 969-340
+49 24 21 969-187
csc@mn-net.com

www.mn-net.com

MACHEREY-NAGEL

MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

Germany and international:
Tel.: +49 24 21 969-0
Fax: +49 24 21 969-199
E-mail: info@mn-net.com

Switzerland:
MACHEREY-NAGEL AG
Tel.: +41 62 388 55 00
Fax: +41 62 388 55 05
E-mail: sales-ch@mn-net.com

France:
MACHEREY-NAGEL EURL
Tel.: +33 388 68 22 68
Fax: +33 388 51 76 88
E-mail: sales-fr@mn-net.com