



SPECIFICATION SHEET

Product Name : Glycine, Extra pure
Product Number : 927.032
Brand : ISOLAB
CAS Number : 56-40-6
Formula : $C_2H_5NO_2$
Formula Weight : 75.06 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	passes test	
Identification	passes test	
Assay	98.5 - 101.5 %	%
Chloride	≤ 0.007 %	%
Sulfate	≤ 0.0065 %	%
Heavy metals	≤ 0.002 %	%
Loss on drying	≤ 0.2 %	%
Residue on ignition	≤ 0.1 %	%
Hydrolyzable substances	passes test	

Quality Management
Isolab Laborgerate GmbH

