



SPECIFICATION SHEET

Product Name : Bismarck brown ®, (C.I. 21010) for microscopy
Product Number : 906.D01
Brand : ISOLAB
CAS Number : 5421-66-9
Formula : $C_{21}H_{24}N_8 \cdot 2HCl$
Formula Weight : 461.39 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Dark brown colored	
Absorption maxima	about 468 nm	nm

Quality Management
Isolab Laborgerate GmbH

