



SPECIFICATION SHEET

Product Name	:	Bismarck brown Y (G), (C.I. 21000) for microscopy
Product Number	:	906.D02
Brand	:	ISOLAB
CAS Number	:	10114-58-6
Formula	:	$C_{18}H_{18}N_8 \cdot 2HCl$
Formula Weight	:	419.31 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Fine black powder	
Absorption maxima	About 457 nm	nm

Quality Management
Isolab Laborgerate GmbH

