



SPECIFICATION SHEET

Product Name : Buffer Solution pH 4.00 ± 0.01 at 20°C
(Potassium Hydrogen Phthalate)
Product Number : 908.B04
Brand : ISOLAB
CAS Number : -
Formula : -
Formula Weight : -

| TEST NAME | SPECIFICATION | UNIT |
|---|---------------------|------|
| pH | 3.99 - 4.01 at 20°C | - |
| Measured with a combination glass electrode after multiple point calibration with reference materials. It is certified traceable to the following National Institute of Standards and Technology (USA), SRM 185i Potassium Hydrogen Phthalate, SRM 186-I-g Potassium Dihydrogen Phosphate and SRM 186-II-g Disodium Hydrogen Phosphate. | | |

Quality Management
Isolab Laborgerate GmbH

