




SPECIFICATION SHEET

Product Name : Buffer Solution pH 7.00 ± 0.01 at 20°C
(Potassium di-hydrogen Phosphate, Disodium Hydrogen Phosphate)
Product Number : 908.B07
Brand : ISOLAB
CAS Number : -
Formula : -
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
pH	6.99 - 7.01 at 20°C	-
Measured with a combination glass electrode after multiple point calibration with reference materials. It is certified traceable to the following National Institute of Standards and Technology (USA), SRM 185i Potassium Hydrogen Phthalate, SRM 186-I-g Potassium Dihydrogen Phosphate and SRM 186-II-g Disodium Hydrogen Phosphate.		





SPECIFICATION SHEET

Product Name : Buffer Solution pH 8.00 ± 0.01 at 20°C
(Di-sodium tetraborate, 1M Hydrochloric Acid)
Product Number : 908.B08
Brand : ISOLAB
CAS Number : -
Formula : -
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
pH	7.99 - 8.01 at 20°C	-
Measured with a combination glass electrode after multiple point calibration with reference materials. It is certified traceable to the following National Institute of Standards and Technology (USA), SRM 191d-I Sodium Bicarbonate, SRM 191d-II Sodium Carbonate, SRM 186-I-g Potassium Dihydrogen Phosphate and SRM 186-II-g Disodium Hydrogen Phosphate.		

Quality Management
Isolab Laborgerate GmbH

