



SPECIFICATION SHEET

Product Name : Buffer Solution pH 8.00 ± 0.01 at 20°C
(Di-sodium tetraborate, 1M Hydrochloric Acid)
Product Number : 908.B08
Brand : ISOLAB
CAS Number : -
Formula : -
Formula Weight : -

| TEST NAME | SPECIFICATION | UNIT |
|---|---------------------|------|
| pH | 7.99 - 8.01 at 20°C | - |
| Measured with a combination glass electrode after multiple point calibration with reference materials. It is certified traceable to the following National Institute of Standards and Technology (USA), SRM 191d-I Sodium Bicarbonate, SRM 191d-II Sodium Carbonate, SRM 186-I-g Potassium Dihydrogen Phosphate and SRM 186-II-g Disodium Hydrogen Phosphate. | | |

Quality Management
Isolab Laborgerate GmbH

