



SPECIFICATION SHEET

Product Name	:	Buffer Solution pH 8.00 ± 0.01 at 20°C (Di-sodium tetraborate, 1M Hydrochloric Acid)
Product Number	:	908.B08
Brand	:	ISOLAB
CAS Number	:	-
Formula	:	-
Formula Weight	:	-

TEST NAME	SPECIFICATION	UNIT
pH	7.99 - 8.01 at 20°C	-
Measured with a combination glass electrode after multiple point calibration with reference materials. It is certified traceable to the following National Institute of Standards and Technology (USA), SRM 191d-I Sodium Bicarbonate, SRM 191d-II Sodium Carbonate, SRM 186-I-g Potassium Dihydrogen Phosphate and SRM 186-II-g Disodium Hydrogen Phosphate.		

Quality Management
Isolab Laborgerate GmbH

