



SPECIFICATION SHEET

Product Name : Buffer Solution pH 9.00 ± 0.01 at 20°C
(Di-sodium Tetraborate, 1M Hydrochloric Acid)
Product Number : 908.B09
Brand : ISOLAB
CAS Number : -
Formula : -
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
pH	8.99 - 9.01 at 20°C	-
Measured with a combination glass electrode after multiple point calibration with reference materials. It is certified traceable to the following National Institute of Standards and Technology (USA), SRM 191d-I Sodium Bicarbonate, SRM 191d-II Sodium Carbonate, SRM 186-I-g Potassium Dihydrogen Phosphate and SRM 186-II-g Disodium Hydrogen Phosphate.		

Quality Management
Isolab Laborgerate GmbH

