

# BRILLIANT GREEN INDICATOR

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 07/02/2017 Version: 0.0


Doc. No: SDS-908.D02/1



### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Product form : Substance  
Trade name : BRILLIANT GREEN INDICATOR  
EC No : 211-190-1  
CAS No : 633-03-4  
Type of product : Pure substance, Hygroscopic substance. Preventive measures apply to the substance in dry state only  
Formula : C<sub>27</sub>H<sub>34</sub>N<sub>2</sub>O<sub>4</sub>S  
Chemical structure :



Synonyms : [4-[4-(diethylamino)benzhydrylene]cyclohexa-2,5-dien-1-ylidene]diethylammonium hydrogen sulphate / 12415 green / ADC brilliant green crystals / aizen diamond green GH / aizen malachite green GH / ammonium, (4-(p-(diethylamino)-alpha-phenylbenzylidene)-2,5-cyclohexadien-1-ylidene) diethyl-, sulfate (1:1) / ammonium, (4-(para-(diethylamino)-alpha-phenylbenzylidene)-2,5-cyclohexadien-1-ylidene) diethyl-, sulfate (1:1) / aniline green / astra diamond green GX / astradiamant green GX / astrazon green D / avon green A-4379 / basic bright green / basic bright green sulfate / basic brilliant green / basic green 1 / basic green I / basic green V / benzaldehyde green / brillant green YNS / brilliant green / brilliant green (hydrogen sulfate) / brilliant green aseptic / brilliant green B / brilliant green BP / brilliant green BP crystals / brilliant green BPC / brilliant green crystals / brilliant green crystals H / brilliant green DSC / brilliant green G / brilliant green GX / brilliant green lake / brilliant green P / brilliant green special / brilliant green sulfate / brilliant green WP crystals / brilliant green Y / brilliant green YN / brilliant green YNS / brilliant lake green Y / brilliant tungstate green toner GT-288 / C.I. 42040 / C.I. basic green 1, sulfate / C.I. basic green 1, sulfate (1:1) / calcozine brilliant green G / deorlene green JJO / diamond green / diamond green G / emerald green / ethanaminium, N-[4-[[4-(diethylamino)phenyl]phenylmethylene]-2,5-cyclohexadien-1-ylidene]-N-ethyl-, sulfate (1:1) / ethyl green (=C.I. basic green 1) / fast green J / fast green JJO / green EN / HIDACO brilliant green / malachite green G / malachite green sulfate / mitsui brilliant green GX / mitsui brilliant green GX / MITSUI brilliant green G / MITSUI brilliant green GX / N-[4-[[4-(diethylamino)phenyl]phenylmethylene]-2,5-cyclohexadien-1-ylidene]-N-ethylethanaminium sulfate (1:1) / solid green / solid green JO / tertrophen brilliant green G / tokyo aniline brilliant green

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

##### 1.2.1. Relevant identified uses

Industrial/Professional use spec : Industrial  
For professional use only  
Use of the substance/mixture : Dye stuff/pigment  
Laboratory chemical

##### 1.2.2. Uses advised against

No additional information available

#### 1.3. Details of the supplier of the safety data sheet

ISOLAB Laborgeräte GmbH  
Am Dillhof 2 - 63863 Eschau / GERMANY  
Tel: + 49 93 74 / 978 55-0  
Fax: +49 93 74 / 978 55-29  
[prodsafe@isolab.de](mailto:prodsafe@isolab.de)

#### 1.4. Emergency telephone number

Country	Organisation/Company	Address	Emergency number	Comment
Germany	Giftnotruf der Charité CBF, Haus VIII (Wirtschaftgebäude), UG	Hindenburgdamm 30 12203 Berlin	+49 30 19240	

# BRILLIANT GREEN INDICATOR

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 2/7/2017

Version: 0.0

Doc. No: SDS-908.D02/1



## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

#### Classification according to Regulation (EC) No. 1272/2008 [CLP]

Acute Tox. 4 (Oral) H302

Eye Irrit. 2 H319

Full text of hazard classes and H-statements : see section 16

#### Adverse physicochemical, human health and environmental effects

No additional information available

### 2.2. Label elements

#### Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :



GHS07

Signal word (CLP) :

Warning

Hazard statements (CLP) :

H302 - Harmful if swallowed  
H319 - Causes serious eye irritation

Precautionary statements (CLP) :

P264 - Wash hands thoroughly after handling  
P280 - Wear eye protection, face protection  
P270 - Do not eat, drink or smoke when using this product.  
P301+P312 - IF SWALLOWED: call a POISON CENTER or doctor/physician if you feel unwell  
P305+P351+P338 - IF IN EYES : Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing  
P337+P313 - If eye irritation persists: Get medical advice/attention  
P501 - Dispose of this material and its container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

### 2.3. Other hazards

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII

This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

## SECTION 3: Composition/information on ingredients

### 3.1. Substances

Name	Product identifier	%
BRILLIANT GREEN INDICATOR	(CAS No) 633-03-4 (EC No) 211-190-1	≈ 95

Full text of H-statements: see section 16

### 3.2. Mixtures

Not applicable

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

First-aid measures general

: Check the vital functions. Unconscious: maintain adequate airway and respiration. Respiratory arrest: artificial respiration or oxygen. Cardiac arrest: perform resuscitation. Victim conscious with laboured breathing: half-seated. Victim in shock: on his back with legs slightly raised. Vomiting: prevent asphyxia/aspiration pneumonia. Prevent cooling by covering the victim (no warming up). Keep watching the victim. Give psychological aid. Keep the victim calm, avoid physical strain. Depending on the victim's condition: doctor/hospital.

First-aid measures after inhalation

: Remove the victim into fresh air. Respiratory problems: consult a doctor/medical service.

First-aid measures after skin contact

: Rinse with water. Soap may be used. Take victim to a doctor if irritation persists.

First-aid measures after eye contact

: Rinse immediately with plenty of water. Take victim to an ophthalmologist if irritation persists.

# BRILLIANT GREEN INDICATOR

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 2/7/2017

Version: 0.0

Doc. No: SDS-908.D02/1



First-aid measures after ingestion : Rinse mouth with water. Immediately after ingestion: give lots of water to drink. Victim is fully conscious: immediately induce vomiting. Call Poison Information Centre ([www.big.be/antigif.htm](http://www.big.be/antigif.htm)). Consult a doctor/medical service if you feel unwell. Ingestion of large quantities: immediately to hospital.

### 4.2. Most important symptoms and effects, both acute and delayed

Symptoms/injuries after inhalation : Slight irritation. Coughing.  
Symptoms/injuries after skin contact : Slight irritation.  
Symptoms/injuries after eye contact : Irritation of the eye tissue.  
Symptoms/injuries after ingestion : Irritation of the gastric/intestinal mucosa. Nausea. Vomiting. Diarrhoea.  
Chronic symptoms : ON CONTINUOUS/REPEATED EXPOSURE/CONTACT: Respiratory difficulties.

### 4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

## SECTION 5: Fire-fighting measures

### 5.1. Extinguishing media

Suitable extinguishing media : Water spray. Polyvalent foam. Alcohol-resistant foam. ABC powder. Carbon dioxide.  
Unsuitable extinguishing media : Solid water jet ineffective as extinguishing medium.

### 5.2. Special hazards arising from the substance or mixture

Fire hazard : DIRECT FIRE HAZARD. Most organic solids may burn if strongly heated.  
Explosion hazard : DIRECT EXPLOSION HAZARD. Most organic solids are liable to dust explosion hazard. INDIRECT EXPLOSION HAZARD. Dust cloud can be ignited by a spark.

### 5.3. Advice for firefighters

Precautionary measures fire : Exposure to fire/heat: keep upwind. Exposure to fire/heat: consider evacuation. Exposure to fire/heat: have neighbourhood close doors and windows.  
Firefighting instructions : Dilute toxic gases with water spray. Take account of environmentally hazardous firefighting water. Use water moderately and if possible collect or contain it.  
Protection during firefighting : Heat/fire exposure: compressed air/oxygen apparatus.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### 6.1.1. For non-emergency personnel

Protective equipment : Gloves. Safety glasses. Protective clothing. Dust cloud production: compressed air/oxygen apparatus.  
Emergency procedures : Mark the danger area. Prevent dust cloud formation. No naked flames. Wash contaminated clothes.

#### 6.1.2. For emergency responders

No additional information available

### 6.2. Environmental precautions

Prevent soil and water pollution. Prevent spreading in sewers.

### 6.3. Methods and material for containment and cleaning up

For containment : Contain released substance, pump into suitable containers. Plug the leak, cut off the supply. Dam up the solid spill. Knock down/dilute dust cloud with water spray.  
Methods for cleaning up : Stop dust cloud by covering with sand/earth. Scoop solid spill into closing containers. Carefully collect the spill/leftovers. Clean contaminated surfaces with a soap solution. Wash clothing and equipment after handling.

### 6.4. Reference to other sections

No additional information available

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

Precautions for safe handling : Comply with the legal requirements. Clean contaminated clothing. Thoroughly clean/dry the installation before use. Do not discharge the waste into the drain. Powdered form: no compressed air for pumping over. Avoid raising dust. Keep away from naked flames/heat. Finely divided: spark- and explosionproof appliances. Finely divided: keep away from ignition sources/sparks. Observe normal hygiene standards. Keep container tightly closed. Measure the concentration in the air regularly. Carry operations in the open/under local exhaust/ventilation or with respiratory protection.

# BRILLIANT GREEN INDICATOR

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 2/7/2017

Version: 0.0

Doc. No: SDS-908.D02/1



### 7.2. Conditions for safe storage, including any incompatibilities

Heat and ignition sources	: KEEP SUBSTANCE AWAY FROM: heat sources. ignition sources.
Information on mixed storage	: KEEP SUBSTANCE AWAY FROM: oxidizing agents. water/moisture.
Storage area	: Store in a dry area. Store at room temperature. Meet the legal requirements.
Special rules on packaging	: SPECIAL REQUIREMENTS: closing. watertight. dry. correctly labelled. meet the legal requirements. Secure fragile packagings in solid containers.
Packaging materials	: SUITABLE MATERIAL: glass.

### 7.3. Specific end use(s)

No additional information available

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

No additional information available

### 8.2. Exposure controls

#### Personal protective equipment:

Protective clothing.

#### Materials for protective clothing:

GIVE EXCELLENT RESISTANCE: No data available. GIVE GOOD RESISTANCE: No data available. GIVE LESS RESISTANCE: No data available. GIVE POOR RESISTANCE: No data available

#### Hand protection:

Gloves

#### Eye protection:

Safety glasses. In case of dust production: protective goggles

Type	Use	Characteristics	Standard
Safety glasses	Dust		

#### Skin and body protection:

Protective clothing

#### Respiratory protection:

Dust production: dust mask with filter type P2

Device	Filter type	Condition	Standard
Dust mask	Type P2	Dust protection	



## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state	: Solid
Appearance	: Crystalline solid. Crystalline powder.
Molecular mass	: 482.64 g/mol
Colour	: Yellow green powder.
Odour	: Odourless.
Odour threshold	: No data available
pH	: No data available
Relative evaporation rate (butylacetate=1)	: No data available
Melting point	: 210 °C
Freezing point	: No data available
Boiling point	: Not applicable

# BRILLIANT GREEN INDICATOR

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 2/7/2017

Version: 0.0

Doc. No: SDS-908.D02/1



Flash point	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: 210 °C
Flammability (solid, gas)	: No data available
Vapour pressure	: No data available
Relative vapour density at 20 °C	: No data available
Relative density	: No data available
Solubility	: Moderately soluble in water. Soluble in ethanol. Water: 1 - 10 g/100ml
Log Pow	: 2.02 (Estimated value)
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: No data available

### 9.2. Other information

VOC content	: 0 %
Other properties	: Translucent. Hygroscopic. Substance has acid reaction.

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

On burning: release of toxic and corrosive gases/vapours (ammonia, nitrous vapours, sulphur oxides, carbon monoxide - carbon dioxide).

### 10.2. Chemical stability

Hygroscopic.

### 10.3. Possibility of hazardous reactions

No additional information available

### 10.4. Conditions to avoid

No additional information available

### 10.5. Incompatible materials

No additional information available

### 10.6. Hazardous decomposition products

No additional information available

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity : Oral: Harmful if swallowed.

BRILLIANT GREEN INDICATOR (633-03-4)	
LD50 oral rat	313 mg/kg (Rat)

Skin corrosion/irritation	: Not classified pH: 2 - 3 (1 %)
Serious eye damage/irritation	: Causes serious eye irritation. pH: 2 - 3 (1 %)
Respiratory or skin sensitisation	: Not classified
Germ cell mutagenicity	: Not classified
Carcinogenicity	: Not classified
Reproductive toxicity	: Not classified
STOT-single exposure	: Not classified
STOT-repeated exposure	: Not classified
Aspiration hazard	: Not classified

# BRILLIANT GREEN INDICATOR

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 2/7/2017

Version: 0.0

Doc. No: SDS-908.D02/1



### SECTION 12: Ecological information

#### 12.1. Toxicity

Ecology - air : Not classified as dangerous for the ozone layer (Regulation (EC) No 1005/2009).  
Ecology - water : Water pollutant (surface water). pH shift. Insufficient data available on ecotoxicity.

#### 12.2. Persistence and degradability

BRILLIANT GREEN INDICATOR (633-03-4)	
Persistence and degradability	Not readily biodegradable in water.

#### 12.3. Bioaccumulative potential

BRILLIANT GREEN INDICATOR (633-03-4)	
Log Pow	2.02 (Estimated value)
Bioaccumulative potential	Low potential for bioaccumulation (Log Kow < 4).

#### 12.4. Mobility in soil

No additional information available

#### 12.5. Results of PBT and vPvB assessment

BRILLIANT GREEN INDICATOR (633-03-4)	
This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII	
This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII	

#### 12.6. Other adverse effects

No additional information available

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

Product/Packaging disposal recommendations : Remove waste in accordance with local and/or national regulations. Hazardous waste shall not be mixed together with other waste. Different types of hazardous waste shall not be mixed together if this may entail a risk of pollution or create problems for the further management of the waste. Hazardous waste shall be managed responsibly. All entities that store, transport or handle hazardous waste shall take the necessary measures to prevent risks of pollution or damage to people or animals. Recycle/reuse. Remove to an authorized incinerator equipped with an afterburner and a flue gas scrubber with energy recovery.

Additional information : LWCA (the Netherlands): KGA category 03. Can be considered as non hazardous waste according to Directive 2008/98/EC.

European List of Waste (LoW) code : 07 03 99 - wastes not otherwise specified

### SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

ADR	IMDG	IATA	ADN	RID
<b>14.1. UN number</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
<b>14.2. UN proper shipping name</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
<b>14.3. Transport hazard class(es)</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
<b>14.4. Packing group</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
<b>14.5. Environmental hazards</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
No supplementary information available				

#### 14.6. Special precautions for user

##### - Overland transport

Not applicable

# BRILLIANT GREEN INDICATOR

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 2/7/2017

Version: 0.0

Doc. No: SDS-908.D02/1



### - Transport by sea

Not applicable

### - Air transport

Not applicable

### - Inland waterway transport

Not applicable

### - Rail transport

Not applicable

### 14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### 15.1.1. EU-Regulations

No REACH Annex XVII restrictions

BRILLIANT GREEN INDICATOR is not on the REACH Candidate List

BRILLIANT GREEN INDICATOR is not on the REACH Annex XIV List

VOC content : 0 %

#### 15.1.2. National regulations

##### Germany

VwVwS Annex reference : Water hazard class (WGK) 2, hazard to waters

WGK remark : Classification water polluting according to external literature source

12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV : Is not subject of the 12. BImSchV (Hazardous Incident Ordinance)

##### Netherlands

SZW-lijst van kankerverwekkende stoffen : The substance is not listed

SZW-lijst van mutagene stoffen : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling : The substance is not listed

##### Denmark

Recommendations Danish Regulation : Young people below the age of 18 years are not allowed to use the product

### 15.2. Chemical safety assessment

No additional information available

## SECTION 16: Other information

Abbreviations and acronyms:

ADN	European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
ATE	Acute Toxicity Estimate
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
PBT	Persistent Bioaccumulative Toxic

# BRILLIANT GREEN INDICATOR

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 2/7/2017

Version: 0.0

Doc. No: SDS-908.D02/1



REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
RID	Regulations concerning the International Carriage of Dangerous Goods by Rail
SDS	Safety Data Sheet
vPvB	Very Persistent and Very Bioaccumulative

Data sources : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

Full text of H- and EUH-statements:	
Acute Tox. 4 (Oral)	Acute toxicity (oral), Category 4
Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
H302	Harmful if swallowed
H319	Causes serious eye irritation

SDS ISOLAB

*This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product*