



SPECIFICATION SHEET

Product Name : Carbol fuchsin powder, for microscopy
Product Number : 909.D01
Brand : ISOLAB
CAS Number : 4197-24-4
Formula : -
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
Appearance	Black powder	
Absorption maximum (%40 Ethanol) λ_{\max}	549-552 nm	nm
Absorptivity (A 1%1 cm in 50% Ethanol, at λ_{\max})	≥ 160	

Quality Management
Isolab Laborgerate GmbH

