



SPECIFICATION SHEET

Product Name : Chloroform, For gas chromatography ECD and FID (Amylene stabilized)
Product Number : 910.054
Brand : ISOLAB
CAS Number : 67-66-3
Formula : CHCl_3
Formula Weight : 119.38 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Clear colorless liquid	
Assay (GC)	$\geq 99.8 \%$	%
Water	$\leq 0.05 \%$	%
Non volatile matter	$\leq 0.0005 \%$	%
Color	$\leq 10 \text{ HU}$	HU
GC/ECD(retention range 1,2,4- trichlorobenzene to decachlorobiphenyl individual signals (lindane standard))	$\leq 3 \text{ pg/ml}$	pg/ml
GC/FID (retention range n-undecane n-tetracontane individual signals (ntetradecane standard))	$\leq 3 \text{ ng/ml}$	ng/ml

Quality Management
Isolab Laborgerate GmbH

