



SPECIFICATION SHEET

Product Name	:	Chloroform, For gas chromatography ECD and FID (Amylene stabilized)
Product Number	:	910.054
Brand	:	ISOLAB
CAS Number	:	67-66-3
Formula	:	CHCl ₃
Formula Weight	:	119.38 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Clear colorless liquid	
Assay (GC)	≥ 99.8 %	%
Water	≤ 0.05 %	%
Non volatile matter	≤ 0.0005 %	%
Color	≤ 10 HU	HU
GC/ECD(retention range 1,2,4- trichlorobenzene to decachlorobiphenyl individual signals (lindane standard))	≤ 3 pg/ml	pg/ml
GC/FID (retention range n-undecane n-tetracontane individual signals (ntetradecane standard))	≤ 3 ng/ml	ng/ml

Quality Management
Isolab Laborgerate GmbH

