



SPECIFICATION SHEET

Product Name : Chloroform, for liquid chromatography , stabilized Amylene
Product Number : 910.057
Brand : ISOLAB
CAS Number : 67-66-3
Formula : CHCl_3
Formula Weight : 119.38 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity (GC)	$\geq 99.8 \%$	%
Identity (IR)	passes test	
Evaporation residue	$\leq 0.0005 \%$	%
Water	$\leq 0.01 \%$	%
Acidity	$\leq 0.0002 \text{ meq/g}$	meq/g
Alkalinity	$\leq 0.0002 \text{ meq/g}$	meq/g
Transmission (at 255 nm)	$\geq 70 \%$	%
Transmission (at 260 nm)	$\geq 85 \%$	%
Transmission (at 300 nm)	$\geq 98 \%$	%
Amylenes as stabilizer	0.003 % - 0.006 %	%

The purpose is only for high pressure liquid chromatography.
Filtered by 0.2 μm filter.

Quality Management
Isolab Laborgerate GmbH

