



SPECIFICATION SHEET

Product Name	:	Chloroform, for liquid chromatography , stabilized Amylene
Product Number	:	910.057
Brand	:	ISOLAB
CAS Number	:	67-66-3
Formula	:	CHCl ₃
Formula Weight	:	119.38 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity (GC)	≥ 99.8 %	%
Identity (IR)	passes test	
Evaporation residue	≤ 0.0005 %	%
Water	≤ 0.01 %	%
Acidity	≤ 0.0002 meq/g	meq/g
Alkalinity	≤ 0.0002 meq/g	meq/g
Transmission (at 255 nm)	≥ 70 %	%
Transmission (at 260 nm)	≥ 85 %	%
Transmission (at 300 nm)	≥ 98 %	%
Amylenes as stabilizer	0.003 % - 0.006 %	%

The purpose is only for high pressure liquid chromatography.
Filtered by 0.2 µm filter.

Quality Management
Isolab Laborgerate GmbH

