



SPECIFICATION SHEET

Product Name : Cyclohexane, Extra pure
Product Number : 912.012
Brand : ISOLAB
CAS Number : 110-82-7
Formula : C₆H₁₂
Formula Weight : 84.16 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity (GC)	≥ 99.8 %	%
Appearance	Colorless, liquid	
Infrared spectrum	passes test	
Residue on evaporation	≤ 0.005 %	%
Non-volatile matter	≤ 0.005%	%

Quality Management
Isolab Laborgerate GmbH

