



### SPECIFICATION SHEET

**Product Name** : N,N-Dimethylformamide, Extra pure  
**Product Number** : 914.042  
**Brand** : ISOLAB  
**CAS Number** : 68-12-2  
**Formula** : C<sub>3</sub>H<sub>7</sub>NO  
**Formula Weight** : 73.09 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity (GC)	≥ 99 %	%
Appearance	Colorless liquid	
Infrared Spectrum	Conforms to structure	
Water	≤ 0.1 %	%
Acidity (HCOOH)	≤ 0.009 %	%

Quality Management  
Isolab Laborgerate GmbH

