



SPECIFICATION SHEET

Product Name : Dibutyl phthalate
Product Number : 914.079
Brand : ISOLAB
CAS Number : 84-74-2
Formula : $C_{16}H_{22}O_4$
Formula Weight : 278.34 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Colourless, liquid	
Infrared spectrum	Conforms to structure	
Specific Gravity at 27 °C	1.047 - 1.049	
Purity (GC)	≥ 99.0 %	%
Color Test (APHA)	≤ 20	
Acidity as Phthalic Acid	≤ 0.01	
Water (by Karl Fischer)	≤ 0.15 %	%

Quality Management
Isolab Laborgerate GmbH

