



SPECIFICATION SHEET

Product Name : Ethylene glycol, Extra pure
Product Number : 920.032
Brand : ISOLAB
CAS Number : 107-21-1
Formula : $C_2H_6O_2$
Formula Weight : 62.07 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity	≥ 99 %	%
Appearance	Colorless, liquid	
Infrared spectrum	passes test	
Color test	≤ 10 APHA	
Acidity	≤ 0.1 ml N %	%
Sulphated Ash	≤ 0.01 %	%

Quality Management
Isolab Laborgerate GmbH

