



SPECIFICATION SHEET

Product Name : Formaldehyde solution about 37%, Extra pure
Product Number : 923.012
Brand : ISOLAB
CAS Number : -
Formula : -
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
Appearance	Clear, colorless	
Assay	36.5 - 38.0 %	%
Methanol	6.0 - 8.0 %	%
Density (25°C)	1.1135-1.095 g/ml	g/ml

Quality Management
Isolab Laborgerate GmbH

