



SPECIFICATION SHEET

Product Name : Giemsa's stain, for microscopy
Product Number : 926.D01
Brand : ISOLAB
CAS Number : 51811-82-6
Formula : -
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
Appearance	Dark green colored powder	
Absorption maxima λ_1 (methanol)	647 - 653 nm	nm
Absorption maxima λ_2 (methanol)	520 - 526 nm	nm
Absorptivity at absorb. max 1	≥ 800	
Absorptivity at absorb. max 2	≥ 350	
Loss on drying (at 105 °C)	$\leq 10 \%$	%

Quality Management
Isolab Laborgerate GmbH

