



SPECIFICATION SHEET

Product Name : n-Hexane, Extra pure
Product Number : 930.022
Brand : ISOLAB
CAS Number : 110-54-3
Formula : C₆H₁₄
Formula Weight : 86.18 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity (GC)	≥ 96.0 %	%
Infrared Spectrum	passes test	
Water	≤ 0.01 %	%
Residue on evaporation	≤ 0.001 %	%
Total sulphur (as S)	≤ 0.01 %	%
Acidity	≤ 0.0005 meq/g	meq/g
Substances darkened by H ₂ SO ₄	passes test	
Cu (Copper)	≤ 0.00002 %	%
Fe (Iron)	≤ 0.00005 %	%
Ni (Nickel)	≤ 0.00002 %	%
Pb (Lead)	≤ 0.00002 %	%

Quality Management
Isolab Laborgerate GmbH

