



### SPECIFICATION SHEET

**Product Name** : Hexadecyltrimethylammonium Bromide, for synthesis  
**Product Number** : 930.05P  
**Brand** : ISOLAB  
**CAS Number** : 57-09-0  
**Formula** :  $C_{16}H_{33}N(CH_3)_3Br$   
**Formula Weight** : 364.45 g/mol

TEST NAME	SPECIFICATION	UNIT
Description	White crystalline powder	
Assay (argentometric)	≥ 97.0 %	%
Identity	passes test	

Quality Management  
Isolab Laborgerate GmbH

