



SPECIFICATION SHEET

Product Name : Hydrogen peroxide 35 %, Extra pure
Product Number : 932.122
Brand : ISOLAB
CAS Number : -
Formula : -
Formula Weight : -

| TEST NAME | SPECIFICATION | UNIT |
|------------------------------------|-------------------------|--------|
| Assay | 34.5 - 36.0 % | % |
| Appearance | Clear colourless liquid | |
| Density (g/cm ³) 25 °C | 1.125 - 1.130 | g/cm3 |
| Color (APHA) | < 10 | |
| Conductivity | > 50 µS/sec | µS/sec |
| Acidity | ≤ 0.03 % | % |
| Stability | ≤ 0.3 ml/min | ml/min |
| Chloride (Cl ⁻) | ≤ 0.0001 % | % |
| Nitrate (NO ₃) | ≤ 0.015 % | % |
| Phosphate (PO ₄) | ≤ 0.005 % | % |
| Sulfate (SO ₄) | ≤ 0.001 % | % |
| As (Arsenic) | ≤ 0.0003 % | % |
| Cu (Copper) | ≤ 0.0001 % | % |
| Fe (Iron) | ≤ 0.00005 % | % |
| Mn (Manganese) | ≤ 0.0001 % | % |
| Pb (Lead) | ≤ 0.001 % | % |
| Sn (Tin) | ≤ 0.001 % | % |
| Non volatile matter | ≤ 0.015 % | % |

Quality Management
Isolab Laborgerate GmbH

