



SPECIFICATION SHEET

Product Name	:	Iodine sublimated, for analysis ACS, ISO, Reag. Ph Eur
Product Number	:	935.016
Brand	:	ISOLAB
CAS Number	:	7553-56-2
Formula	:	I ₂
Formula Weight	:	253.8 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay	≥ 99.8 %	%
Identity	passes test	
Bromide and chlorine (as Cl)	≤ 0.005 %	%
Residue	≤ 0.01 %	%

Quality Management
Isolab Laborgerate GmbH

