



SPECIFICATION SHEET

Product Name : May and Grunwald's stain, for microscopy
Product Number : 946.D02
Brand : ISOLAB
CAS Number : -
Formula : -
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
Appearance	A dark green powder	
Loss on drying (at 105°C)	≤ 10 %	%
Absorbance maxima at λ1	648-653 nm	nm
Absorbance maxima at λ2	520 - 525 nm	nm
Specific absorptivity A 1%/1 cm at λ1	1100 - 1300	
Specific absorptivity A 1%/1 cm at λ2	600 - 800	

Quality Management
Isolab Laborgerate GmbH

