



SPECIFICATION SHEET

Product Name : Nigrosine alcohol soluble, (C.I. 50415) for microscopy
Product Number : 951.D02
Brand : ISOLAB
CAS Number : 11099-03-9
Formula : -
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
Appearance	Black powder	
Absorption maxima	about 560 nm	nm
Specific absorptivity (A1%,1 cm 0.025 g/L ethanol 50%)	200 - 300	

Quality Management
Isolab Laborgerate GmbH

