



SPECIFICATION SHEET

**Product Name** : Oil red O, (C.I. 26125) for electrophoresis  
**Product Number** : 955.D01  
**Brand** : ISOLAB  
**CAS Number** : 1320-06-5  
**Formula** :  $C_{26}H_{24}N_4O$   
**Formula Weight** : 408.5 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Brownish red powder	
Absorption maxima( $\lambda$ ) (Chloroform)	528 - 534 nm	nm
Loss on drying 105°C	≤ 5 %	%

Quality Management  
Isolab Laborgerate GmbH

