



SPECIFICATION SHEET

Product Name : Orcein, For microscopy
Product Number : 957.D02
Brand : ISOLAB
CAS Number : 1400-62-0
Formula : -
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
Appearance	Dark brown to black-purplish powder	
Absorption in NaOH 0.01 mol/L	575 - 582 nm	nm
Specific absorptivity	≥ 500	
Loss on drying	≤ 7 %	%

Quality Management
Isolab Laborgerate GmbH

