



### SPECIFICATION SHEET

**Product Name** : n-Pentane, Extra pure  
**Product Number** : 959.012  
**Brand** : ISOLAB  
**CAS Number** : 109-66-0  
**Formula** : C<sub>5</sub>H<sub>12</sub>  
**Formula Weight** : 72.15 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity (GC)	≥ 99.0 %	%
Identity (IR)	passes test	
Free acid (as CH <sub>3</sub> COOH)	≤ 0.001 %	%
Non volatile matter	≤ 0.005 %	%
Readily carbonizable substances	passes test	

Quality Management  
Isolab Laborgerate GmbH

