



SPECIFICATION SHEET

Product Name : Petroleum Ether boiling range 40-60°C, DAB
Product Number : 959.033
Brand : ISOLAB
CAS Number : -
Formula : -
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
Appearance	Clear colorless liquid	
Acidity or alkalinity	passes test	
Boiling range (40-60°C)	≥ 75.0 % (v/v)	%
Colour (APHA)	≤ 10	
Foreign odour, non-volatile matter	passes test	
Sulphur compounds, reducing substances	passes test	
Benzene (GC)	≤ 0.001 %	%
n-Hexane (GC)	≤ 2.0 %	%
Tetraethyllead	passes test	
Readily carbonisable substance	passes test	
Evaporation residue	≤ 0.001 %	%
Sulphur compounds (as S)	≤ 0.005 %	%
Substances dark by H ₂ SO ₄	passes test	
Water	≤ 0.01 %	%

Quality Management
Isolab Laborgerate GmbH

