



SPECIFICATION SHEET

Product Name : Potassium chromate, For analysis ACS, Reag. Ph Eur
Product Number : 960.046
Brand : ISOLAB
CAS Number : 7789-00-6
Formula : CrK_2O_4
Formula Weight : 194.19 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay	99.5 - 100.5 %	%
In water insoluble matter	≤ 0.003 %	%
pH-value (5%; water)	9.0 - 9.8	
Chloride (Cl)	≤ 0.001 %	%
Sulphate (SO_4)	≤ 0.01 %	%
Ca (Calcium)	≤ 0.005 %	%
Na (Sodium)	≤ 0.02 %	%
Pb (Lead)	≤ 0.005 %	%
Loss on drying	≤ 0.5 %	%

Quality Management
Isolab Laborgerate GmbH

