



SPECIFICATION SHEET

Product Name : Potassium dichromate, For analysis ACS, ISO, Reag. Ph Eur
Product Number : 960.056
Brand : ISOLAB
CAS Number : 7778-50-9
Formula : $\text{Cr}_2\text{K}_2\text{O}_7$
Formula Weight : 294.19 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (Iodometric, calc. on dried substance)	$\geq 99.9 \%$	%
Insoluble matter	$\leq 0.003 \%$	%
Chloride (Cl)	$\leq 0.001 \%$	%
Sulphate (SO_4)	$\leq 0.005 \%$	%
Ca (Calcium)	$\leq 0.002 \%$	%
Cu (Copper)	$\leq 0.001 \%$	%
Fe (Iron)	$\leq 0.001 \%$	%
Na (Sodium)	$\leq 0.02 \%$	%
Pb (Lead)	$\leq 0.005 \%$	%
Loss on drying (105 °C)	$\leq 0.05 \%$	%

Quality Management
Isolab Laborgerate GmbH

