



SPECIFICATION SHEET

Product Name : Potassium hydroxide, Pellets pure
Product Number : 960.082
Brand : ISOLAB
CAS Number : 1310-58-3
Formula : KOH
Formula Weight : 56.11 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (KOH)	≥ 85.0 %	%
Carbonate (as K ₂ CO ₃)	≤ 2.0 %	%
Chloride (Cl)	≤ 0.01 %	%
Sulphate (SO ₄)	≤ 0.01 %	%
Heavy metals (as Pb)	≤ 0.002 %	%
Al (Aluminium)	≤ 0.002 %	%
Fe (Iron)	≤ 0.002 %	%

Quality Management
Isolab Laborgerate GmbH

