



### SPECIFICATION SHEET

**Product Name** : Polysorbate 80, For synthesis  
**Product Number** : 960.161  
**Brand** : ISOLAB  
**CAS Number** : 9005-65-6  
**Formula** : -  
**Formula Weight** : -

| TEST NAME             | SPECIFICATION | UNIT |
|-----------------------|---------------|------|
| Density (d 20°C/ 4°C) | 1.060 - 1.090 |      |
| Saponification value  | 45 - 55       |      |
| Hydroxyl value        | 65 - 80       |      |
| Identity (IR)         | passes test   |      |

Due to its specific melting range the product may be solid, liquid, a solidified melt or a supercooled melt.

Quality Management  
Isolab Laborgerate GmbH

