



### SPECIFICATION SHEET

**Product Name** : Potassium Permanganate, Extra pure  
**Product Number** : 960.172  
**Brand** : ISOLAB  
**CAS Number** : 7722-64-7  
**Formula** : -  
**Formula Weight** : -

TEST NAME	SPECIFICATION	UNIT
Assay (manganometric)	99.0 - 100.5 %	%
Identity	passes test	
Appearance of solution	passes test	
In water insoluble matter	≤ 0.2 %	%
Chloride (Cl)	≤ 0.01 %	%
Sulfate (SO <sub>4</sub> )	≤ 0.03 %	%

Quality Management  
Isolab Laborgerate GmbH

