



PRODUCT SPECIFICATION

Product Name : 2-Propanol $\geq 99.8\%$ Extra Pure
Product Number : 961.022
Brand : ISOLAB
CAS Number : 67-63-0
Formula : C_3H_8O
Formula Weight : 60.1 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity	≥ 99.8	%
Water	≤ 0.1	%
Acidity as acetic acid	≤ 0.001	%
Water Miscibility	Miscible	-
Nonvolatile Residue	≤ 0.001	g/100ml
Appearance	Clear & Free From Suspended Matter	-
Color, Pt-Co	Max 5.0	Pt Co
Density @20C	0.785 - 0.786	g/ml
Distillation, IBP	Min. 81.8	$^{\circ}C$
Distillation, DP	Max. 82.8	$^{\circ}C$
Refractive Index @20C	1.376 - 1.378	-

Quality Management
Isolab Laborgerate GmbH

