



SPECIFICATION SHEET

Product Name : 2-Propanol, For gas chromatography ECD and FID
Product Number : 961.024
Brand : ISOLAB
CAS Number : 67-63-0
Formula : C₃H₈O
Formula Weight : 60.1 g/mol

| TEST NAME | SPECIFICATION | UNIT |
|--|------------------------|-------|
| Appearance | Clear colorless liquid | |
| Assay (GC) | ≥ 99.8 % | % |
| Color | ≤ 10 HU | HU |
| Water (KF) | ≤ 0.1 % | % |
| Non volatile matter | ≤ 0.0005 % | % |
| GC/ECD (retention range 1,2,4-trichlorobenzene to decachlorobiphenyl individual signal (lindane standard)) | ≤ 3 pg/ml | pg/ml |
| GC/FID (retention range n-undecane - n-tetracontane individual signals (n-tetradecane standard)) | ≤ 3 ng/ml | ng/ml |

Quality Management
Isolab Laborgerate GmbH

