

SPECIFICATION SHEET

Product Name : Pyronin B, (C.I.45010)
Product Number : 961.D01
Brand : ISOLAB
CAS Number : 2150-48-3
Formula : $C_{42}H_{54}Cl_8Fe_2N_4O_2$
Formula Weight : 1042.22 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Dark green powder	
Absorption maxima	about 553 nm	nm
Specific Absorptivity A1% 1 cm (λ_{max} ; 0.0003%, ethanol 50%)	≥ 1800	

Quality Management
Isolab Laborgerate GmbH


