



SPECIFICATION SHEET

Product Name : Salicylic acid, Extra pure
Product Number : 967.012
Brand : ISOLAB
CAS Number : 69-72-7
Formula : $C_7H_6O_3$
Formula Weight : 138.12 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	White, crystalline, powder	
Assay	$\geq 99.5 \%$	%
Residue on ignition	$\leq 0.1 \%$	%
Heavy metals	$\leq 0.002 \%$	%
Chloride (Cl)	$\leq 0.01 \%$	%
Sulfate	$\leq 0.03 \%$	%
Loss on drying	$\leq 0.5 \%$	%
Phenol	$\leq 0.05 \%$	%

Quality Management
Isolab Laborgerate GmbH

