



SPECIFICATION SHEET

Product Name : Silver sulfate extra pure
Product Number : 968.052
Brand : ISOLAB
CAS Number : 10294-26-5
Formula : Ag_2SO_4
Formula Weight : 311.8 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (argentometric)	≥ 98.5 %	%
Identity	passes test	
Chloride (Cl)	≤ 0.002 %	%
Cu (Copper)	≤ 0.005 %	%
Fe (Iron)	≤ 0.005 %	%
Pb (Lead)	≤ 0.005 %	%
Substances not precipitated by hydrochloric Acid (as SO_4)	≤ 0.03 %	%

Quality Management
Isolab Laborgerate GmbH

