



SPECIFICATION SHEET

Product Name : Silver Sulfate GR for analysis ACS
Product Number : 968.055
Brand : ISOLAB
CAS Number : 10294-26-5
Formula : Ag_2SO_4
Formula Weight : 311.8 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (argentometric)	$\geq 98.5 \%$	%
Identity	passes test	
Insoluble matter and silver chloride	$\leq 0.02 \%$	%
Chloride (Cl)	$\leq 0.001 \%$	%
Nitrate (NO_3)	$\leq 0.001 \%$	%
Cu (Copper)	$\leq 0.0005 \%$	%
Fe (Iron)	$\leq 0.001 \%$	%
Pb (Lead)	$\leq 0.001 \%$	%
Substances not precipitated by hydrochloric Acid (as SO_4)	$\leq 0.03 \%$	%

Quality Management
Isolab Laborgerate GmbH

