



### SPECIFICATION SHEET

Product Name	:	Sodium Disulfite Ph Eur,BP,JP,NF,E 223
Product Number	:	969.073
Brand	:	ISOLAB
CAS Number	:	7681-57-4
Formula	:	Na <sub>2</sub> O <sub>5</sub> S <sub>2</sub>
Formula Weight	:	190.11 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub> )	≥ 98.0 - 99.5	%
Assay (SO <sub>2</sub> )	≥ 66.0 - 67.0	%
Identity	passes test	
Appearance of solution	passes test	
pH (10% sol)	4.0 - 4.5	
Chloride (Cl)	≤ 0.005	%
Heavy metals (as Pb)	≤ 0.001	%
As (Arsenic)	≤ 0.00005	%
Fe (Iron)	≤ 0.0005	%
Hg (Mercury)	≤ 0.00001	%
Pb (Lead)	≤ 0.00005	%
Se (Selenium)	≤ 0.00005	%

Quality Management  
Isolab Laborgerate GmbH

