



SPECIFICATION SHEET

Product Name : Sodium sulfite anhydrous, Ph Eur, BP, NF, E 221
Product Number : 969.153
Brand : ISOLAB
CAS Number : 7757-83-7
Formula : Na_2SO_3
Formula Weight : 126.04 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (Na_2SO_3)	$\geq 98.0\%$	%
Assay (SO_2)	$\geq 49.8\%$	%
Identity	passes test	
Appearance of solution	passes test	
Thiosulfate (S_2O_3)	$\leq 0.1\%$	%
Heavy metals (as Pb)	$\leq 0.001\%$	%
As (Arsenic)	$\leq 0.0001\%$	%
Hg (Mercury)	$\leq 0.00001\%$	%
Fe (Iron)	$\leq 0.001\%$	%
Pb (Lead)	$\leq 0.0001\%$	%
Se (Selenium)	$\leq 0.0001\%$	%
Zn (Zinc)	$\leq 0.0001\%$	%

Quality Management
Isolab Laborgerate GmbH

