

SPECIFICATION SHEET

Product Name : Sodium thiosulfate anhydrous
Product Number : 969.162
Brand : ISOLAB
CAS Number : 7772-98-7
Formula : $\text{Na}_2\text{O}_3\text{S}_2$
Formula Weight : 158.11 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (Iodometric, on dried substance)	≥ 98.0 %	%
Identity	passes test	
pH-value (5%, water)	6.0 - 8.5	
Insoluble matter	≤ 0.01 %	%
Sulphate and sulphite (SO_4)	≤ 0.5 %	%
Sulphide (S)	≤ 0.001 %	%
Calcium (Ca)	≤ 0.005 %	%
Heavy metals (as Pb)	≤ 0.002 %	%
Fe (Iron)	≤ 0.002 %	%
Loss on drying	≤ 1.0 %	%

Quality Management
Isolab Laborgerate GmbH


