



### SPECIFICATION SHEET

Product Name : Starch, Soluble extra pure  
Product Number : 970.012  
Brand : ISOLAB  
CAS Number : 9005-84-9  
Formula :  $(C_6H_{10}O_5)_n$   
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
Appearance	Almost white powder	
Sensitivity with Iodine solution	passes test	
Reducing substances	passes test	
Microbiological test	passes test	
Colony count (aerobic bacteria)	$\leq 5000$ CFU/g	CFU/g
Colony count (Yeasts and moulds)	$\leq 500$ CFU/g	CFU/g
Salmonella species	absent in 10 g	g
Staphylococcus aureus	absent in 1 g	g
Pseudomonas aeruginosa	absent in 1 g	g
E.coli	absent in 1 g	g

Quality Management  
Isolab Laborgerate GmbH

