



SPECIFICATION SHEET

Product Name : Toluene $\geq 99.5\%$, for analysis ACS
Product Number : 973.015
Brand : ISOLAB
CAS Number : 108-88-3
Formula : C_7H_8
Formula Weight : 92.14 g/mol

| TEST NAME | SPECIFICATION | UNIT |
|-------------------------------------|-----------------------|------|
| Appearance (Color) | Colorless | |
| Appearance (Form) | Liquid | |
| Color Test | < 10 APHA | APHA |
| Sulfur Compounds (as S) | ≤ 0.003 % | % |
| Purity (GC) | ≥ 99.50 % | % |
| Water (by Karl Fischer) | ≤ 0.03 % | % |
| Infrared spectrum | Conforms to Structure | |
| Residue on Evaporation | ≤ 0.001 % | % |
| Substances Darkened (by H_2SO_4) | Passes test | |
| Benzene vol % | ≤ 0.005 % | % |
| Acidity mg NaOH/100ml | Not detected | |

Quality Management
Isolab Laborgerate GmbH

