



SPECIFICATION SHEET

Product Name : Trypan blue, (C.I.23850) for microscopy
Product Number : 974.D01
Brand : ISOLAB
CAS Number : 72-57-1
Formula : $C_{34}H_{24}N_6Na_4O_{14}S_4$
Formula Weight : 960.79 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Dark green to brown or brown-black powder, Crystalline powder or granule	
Dye content (spectrophotometrically)	≥ 70 %	%

Quality Management
Isolab Laborgerate GmbH

