



SPECIFICATION SHEET

Product Name : Victoria blue, (C.I.44045) for microscopy
Product Number : 980.D01
Brand : ISOLAB
CAS Number : 2580-56-5
Formula : $C_{33}H_{32}ClN_3$
Formula Weight : 506.08 g/mol

| TEST NAME | SPECIFICATION | UNIT |
|-----------------|---------------------|------|
| Appearance | Brownish red powder | |
| Absorption max. | About 600 nm | nm |

Quality Management
Isolab Laborgerate GmbH

