



### SPECIFICATION SHEET

**Product Name** : Wright's stain, For microscopical staining  
**Product Number** : 986.D01  
**Brand** : ISOLAB  
**CAS Number** : 68988-92-1  
**Formula** : -  
**Formula Weight** : -

TEST NAME	SPECIFICATION	UNIT
Appearance	Greenish black/blue crystalline powder	
Absorption max.	About 600 nm	nm
Absorption maxima in methanol ( $\lambda_1$ )	640 - 648 nm	nm
Absorption maxima in methanol ( $\lambda_2$ )	522- 525 nm	nm
Specific absorptivity at $\lambda_1$	$\geq 1100$	
Specific absorptivity at $\lambda_2$	$\geq 650$	
Loss on drying (at 110°C)	$\leq 10 \%$	%
Suitability for microscopy	passes test	

Quality Management  
Isolab Laborgerate GmbH

